

| | | | |
|--|------------------------------------|---------------------------------|--------|
| Document Title: Slick, AquaLab Pre | | Part # and Rev. 18078 | |
| | | Release Date: | |
| Rev. | Description | Revision By | Date |
| Date -- Time | Added description from other slick | Allison | 6/8/16 |

Production Filename:

<http://publications.decagon.com/Marketing/Description%20Files/18078 Slick AquaLab Pre.pdf>

Printing Process: Digital

Size: 11 inches wide, 8.5 inches tall

Ink: 4/4 Color Both Sides

Paper: 100 lbs. book gloss coated

Page Count: 2

Proof type: soft proof

Single hole drilling: none

WATER ACTIVITY METER



AquaLab PRE

- **Accurate** —
±0.01 a_w
- **Verifiable** —
with independent salt standards.
- **Repeatable** —
different users, different locations, same result.
- **Easy to use** —
precise measurements with minimal training.

AFFORDABLE BENCHTOP WATER ACTIVITY

Repeatable measurements, solid support

The AquaLab Pre is a robust entry-level water activity meter with stripped down form, not limited function. It uses the same industry-standard primary a_w method you find in the top-of-the-line AquaLab Series 4TE. The Pre may look basic, but that's because all the good stuff is inside.

Minimum maintenance. The Pre uses the dewpoint method—a primary measurement of water

activity, that eliminates the need for lengthy sensory calibration. **Just a quick verification check** and you're good to go.

Real water activity—FAST. No need to sacrifice accuracy for speed. Primary dewpoint measurement allows the Pre to deliver **top-accuracy readings in 5 minutes** or less.

Solid repeatability. Water activity is temperature

dependent, measuring at the same temperature every time assures consistency in your readings. The Pre holds the sample at 25 °C so **temperature fluctuations won't affect your readings.**

Reliable accuracy. The Pre dewpoint method and temperature-stable features mean good agreement between the Pre and the Series 4. Its basic **0.01 a_w accuracy** enables many of our customers to add water

activity testing at the line, loading dock, or offsite facility.

Dependable support. The Pre includes Deacon's **free technical support**—food scientists and instrument specialists just a phone call away.

Full range. The Pre lets you measure over the whole range of water activities, **0.05 to 1.00 a_w .**

AQUA[®] LAB
BY DECA[®]GON

2365 NE Hopkins Court
Pulman, Washington 99163
509-332-7756
sales@aqualab.com
www.aqualab.com