

Document Title: Slick, Hyprop		Part # and Rev. 18091	
		Release Date:	
Rev.	Description	Revision By	Date
Date -- Time	Added description from Sharpdots	Allison	6/9/16

Production Filename:

http://publications.decagon.com/Marketing/Description%20Files/18091_Slick_Hyprop.pdf

Printing Process: Digital

Size: 8.5 inches wide, 11 inches wide

Ink: 4/4 Color Front Only

Paper: 100 lbs Book Gloss Coated

Page count: 2

Proof type: soft proof

Single Hole Drilling: none



HYPROP



Generate Moisture Release Curves

Forget the pain of using pressure plates. Use Hyprop to generate detailed moisture release curves, typically within a day.

Just take your undisturbed sample, insert the measurement head, put it on the balance, and one to two days, you have a detailed moisture release curve.

- *Uses undisturbed soil samples.*
- *Works in all soil types*
- *Bonus feature gives unsaturated hydraulic conductivity values for the soil sample.*
- *Range: 0 to -85 kPa*