

Document Title: Slick, Pawkit		Part # and Rev. 18102	
		Release Date:	
Rev.	Description	Revision By	Date
Date -- Time	Added description from other slicks	Allison	6/8/16

Production Filename:

http://publications.decagon.com/Marketing/Description%20Files/18102_Slick_Pawkit.pdf

Printing Process: Digital

Size: 8.5 inches wide, 11 inches tall

Ink: 4/4 Color Both Sides

Paper: 100 lbs. book gloss coated

Page Count: 2

Proof type: soft proof

Single hole drilling: none



Ultra-compact water activity kit.

- Low-maintenance sensor and readout in a single unit.
- Perform quick checks, measure a_w in 5 minutes or less.
- Easy-to-read display shows water activity and sample temperature together.
- Long 3-year battery life. No adapter or external power supply required.
- Remarkably portable, and lightweight.
- One-year warranty—parts and labor. Extended warranty available.
- Easy calibration saves time.
- Included: protective case, caps, solvents, quick-start guide, and manual.

■ Designed for quick checks by health inspectors, food scientists, and QC managers in foods, pharmaceuticals, and cosmetics.

Water activity—four inches sleek and four ounces slim.

When you can't bring the lab with you. The new portable Pawkit water activity kit is a lightweight, space-saving device for any on-the-go Quality Control Professional.

Four inches. This sleek device is compact, lightweight, portable and discreet. At four inches long, it takes up no more room than a small calculator.

Four ounces. The ultra-compact Pawkit delivers $\pm 0.02 a_w$ measurements in 5 minutes, perfect for quick checks of facilities without laboratories. ■

Pawkit
AquaLab

AQUA LAB
BY DESIGN

2365 NE Hoggins Court
Fullman, Washington 99163
509-332-2756
sales@aqualab.com
www.aqualab.com