

Document Title: <b>Pawkit Upgrade Kit Insert</b>		Part # <b>13774</b>	
		Release Date: <b>06/01/11</b>	
Rev.	Description	Revision By	Date
Date -- Time	Uploaded insert to repository. Available at <a href="http://publications.decagon.com">http://publications.decagon.com</a> . Please ask archivist for previous versions or use Beanstalk application.	NJR	10/29/13
	Changed measurement order from height x width to width x height	Allison	5/23/2016

**Decagon and Vendor Production File Name:**

[http://manuals.decagon.com/Inserts/13774\\_Pawkit Upgrade Kit\\_Print.pdf](http://manuals.decagon.com/Inserts/13774_Pawkit Upgrade Kit_Print.pdf)

**Dimensions:** 5 inches wide, 2.75 inches tall

**Colors:** Full Color

**Material:** 80 lbs/ 216 g/m<sup>2</sup> Long Grain White (14583)

**Finish:** Cut to Size

**Adhesive:** None

**Serialization:** None

**Serial Number Example:** NA

**Serial Number Text:** NA

**Special Notes:** Pages are duplexed front to back

Illustration is Ref Only

**\*\*Not to scale \*\***



*Pawkit Upgrade Kit*

This upgrade kit contains all parts needed to upgrade an AquaLab Pawkit from a GE Panametrics MiniCap 2 (MC-2) sensor to a Honeywell (HCH-1000-001) RH sensor. The new sensor and PCBA in this kit have been fully tested for functionality but are uncalibrated.

Follow the battery replacement instructions to exchange out old parts with the new parts included in this kit. After re-assembly run a full Pawkit calibration.

*Upgrade Kit Is For Authorized Distributors Only*