

Document Title: Description, Label, AquaLab 4TE		Part # 13978	
		Release Date: 8/31/2010	
Rev.	Description	Revision By	Date
01	File Created	Tyler J Zollinger	8/31/2010

Production Filename: 13978 (In Product Library)

File Path: Decadoc\AquaLab Series 4\Marketing\13978 Label, AquaLab 4TE

Dimensions: 4.25 inch wide, 1.25 inch tall with corner radius of .0938

Material: Urethane domed label

Colors: Pantone Reflex Blue, Pantone PMS 427C, Pantone 320

Adhesive: Permanent acrylic adhesive

Special Notes: Illustrations are Ref Only ** Not to Scale **

